





# Mindray Neutralizing Antibody Assays Information



# *Neutralizing antibodies and Vaccines*




# Neutralizing antibodies bind to the protein, and prevent recognition of the ACE2 receptor and or viral fusion



*SARS-CoV-2 Neutralizing Antibody Mechanism*

Image from:  
<https://www.prosci-inc.com/ace2-antibodies/>

# Leading vaccine candidates against SARS-COV-2



|   |  |  |   |   |
|---|--|--|---|---|
| <b>COVID-19 vaccines using this platform</b>                            | Sinopharm  Sinovac<br>Bharat Biotech and etc.  | CanSino  Oxford-AstraZeneca<br>Gamaleya  Johnson & Johnson   | Pfizer-BioNTech  Moderna<br>and etc.  | Vector Institute  Novavax<br>Anhui Zhifei Longcom and etc.  |
| <b>Phase and status</b>   | <ul style="list-style-type: none"> <li>● <b>Sinopharm-Phase 3</b>-Approved in China, U.A.E., Bahrain. Emergency use in other countries.</li> <li>● <b>Sinovac-Phase 3</b>-Approved in China. Emergency use in other countries.</li> <li>● <b>Bharat Biotech-Phase 3</b>-Emergency use in India, Iran, Zimbabwe.</li> </ul> | <ul style="list-style-type: none"> <li>● <b>CanSino-Phase 3</b>-Approved in China. Emergency use in Mexico, Pakistan.</li> <li>● <b>AstraZeneca-Phase 2/3</b>-Approved in Brazil. Emergency use in U.K., E.U., other countries.</li> <li>● <b>Gamaleya-Phase 3</b>-Early use in Russia. Emergency use in other countries.</li> <li>● <b>J&amp;J-Phase 3</b>-Emergency use in U.S., E.U., Bahrain.</li> </ul> | <ul style="list-style-type: none"> <li>● <b>Pfizer-BioNTech-Phase 2/3</b>-Approved in several countries. Emergency use in U.S., E.U., other countries.</li> <li>● <b>Moderna-Phase 3</b>-Approved in Switzerland. Emergency use in U.S., U.K., E.U., others.</li> </ul> | <ul style="list-style-type: none"> <li>● <b>Vector Institute-Phase 3</b>-Early use in Russia.</li> <li>● <b>Novavax-Phase 3</b>.</li> <li>● <b>Anhui Zhifei Longcom-Chinese Academy of Medical Sciences-Phase 3</b>-Emergency use in China and Uzbekistan.</li> </ul> |
| <b>Mindray' solution for monitoring vaccine-induced immune response</b> | <i>Anti-S-RBD IgG/Total antibody</i>   | <i>Neutralizing Antibody/Anti-S-RBD IgG/Total antibody??</i>   | <i>Neutralizing Antibody/Anti-S-RBD IgG/Total antibody??</i>  | <i>Anti-S-RBD IgG/Total antibody??</i>  |

# How the inactivated vaccine Works

| Vaccine developer | Number of doses | Dosing schedule(day) | Trial group    | Monitoring time of neutralizing antibody(day)  | Efficacy                                  |
|-------------------|-----------------|----------------------|----------------|--|---|
| Sinopharm         | 2               | 0+21 (Phase 2)       | 1. 18–59 years | PRNT:0+35+49 (peak-28 days after the second inoculation) <sup>1</sup><br>Neutralizing antibody GMT: 282.7149<br>specific IgG-binding antibody:0+35 (peak-14 days after the second inoculation) <sup>2</sup><br>IgG titers (to whole SARS-CoV-2 antigen): 215   | 79.34%                                    |
|                   |                 | 0+28 (Phase 1)       | 2. ≥60 years   | PRNT:0+7+14+28+32+42 (peak-14 days after the second inoculation) <sup>1</sup><br>Neutralizing antibody GMT: 131.4649; Seroconversion rate: 100%  |   |
| Sinovac Biotech   | 2               | 0+14 (Phase 2)       | 1. 18–59 years | PRNT:0+28+42 (peak-14 days after the second inoculation) <sup>3</sup><br>Neutralizing antibody GMT: 27.6; Seroconversion rate: 92.4%<br>RBD specific IgG:0+28+42 (peak-14 days after the second inoculation)<br>IgG titers: 1094.3; Seroconversion rate: 96.5% | 50.38%-Brazil trial<br>83.5%-Turkey trial |
|                   |                 | 0+28 (Phase 1)       | 2. ≥60 years   | PRNT:0+28+56 (peak-28 days after the second inoculation) <sup>4</sup><br>Neutralizing antibody GMT: 54.9; Seroconversion rate: 100.0%  |   |

# How the mRNA vaccine Works

| Vaccine developer | Number of doses | Dosing schedule(day) | Trial group    | Monitoring time of neutralizing antibody(day)   | Human convalescent serum   | Efficacy |
|-------------------|-----------------|----------------------|----------------|---|--|----------|
| Moderna           | 2               | 0+28                 | 1. 18-55 years | Live virus PRNT <sub>80</sub> :0+42 (peak-14 days after the second inoculation) <sup>5</sup> GMR: 654.3<br>IgG ELISA-S-2P:0+14+28+35+42+56 (peak-14 days after the second inoculation) GMT: 811,119<br>IgG ELISA-RBD:0+14+28+35+42+56 (peak-14 days after the second inoculation) GMT: 558,905          | Live virus PRNT <sub>80</sub> :158.3<br>IgG ELISA-S-2P:142,140<br>IgG ELISA-RBD:37,857 | 94.10%   |
|                   |                 |                      | 2. 56-70 years | Live virus PRNT <sub>80</sub> :0+42 (peak-14 days after the second inoculation) <sup>6</sup> GMR: 878<br>IgG ELISA-S-2P:0+14+28+35+42+56 (peak-7 days after the second inoculation) GMT: 5,033,017<br>IgG ELISA-RBD:0+14+28+35+42+56 (peak-7 days after the second inoculation) GMT: 1,471,882          |  |          |
|                   |                 |                      | 3. ≥71 years   | Live virus PRNT <sub>80</sub> :0+42 (peak-14 days after the second inoculation) <sup>6</sup> GMR: 317<br>IgG ELISA-S-2P:0+14+28+35+42+56 (peak-7 days after the second inoculation) GMT: 2,636,979<br>IgG ELISA-RBD:0+14+28+35+42+56 (peak-7 days after the second inoculation) GMT: 711,752            |  |          |
| Pfizer-BioNTech   | 2               | 0+21                 | 1. 18-55 years | PRNT:0+7+21+28+42+49+84 (peak-7 days after the second inoculation) <sup>7</sup> Neutralization Titer:312<br>S1 specific IgG:0+7+21+28+42+49+84 (peak-7 days after the second inoculation) (U/ml): 8,279<br>RBD specific IgG:0+7+21+28+42+49+84 (peak-7 days after the second inoculation) (U/ml): 6,299 | PRNT50-Neutralization Titer:94<br>S1-IgG (U/ml): 631<br>RBD-IgG (U/ml):602             | 91.30%   |
|                   |                 |                      | 2. 65-85 years | PRNT <sub>50</sub> :0+21+28+35 (peak-14 days after the second inoculation) <sup>8</sup> Neutralization Titer: 206<br>S1-Binding IgG:0+21+28+35 (peak-7 days after the second inoculation) (U/ml): 7,985   |  |          |

# How the Viral vector-based vaccine Works

| Vaccine developer            | Number of doses | Dosing schedule(day) | Trial group    | Monitoring time of neutralizing antibody(day)  | Human convalescent serum   | Efficacy   |
|------------------------------|-----------------|----------------------|----------------|--|--|--|
| <b>Gamaleya</b>              | 2               | 0+21                 | 18–60 years    | PRNT:0+14+28+42 (peak-21 days after the second inoculation) <sup>9</sup><br>GMT: 49.25; Seroconversion rate:100%<br>RBD specific IgG:0+14+21+28+42 (peak-21 days after the second inoculation)<br>GMT: 14703; Seroconversion rate:100%   | at ~1 month after recovery<br>PRNT- GMT: 32.96;<br>Seroconversion rate:90.9%<br>RBD IgG- GMT: 1266;<br>Seroconversion rate:85.8% | 91.60%   |
| <b>AstraZeneca</b>           | 2               | 0+28                 | 1. 18-55 years | live SARS-CoV-2 micro-neutralization assay MNA <sub>80</sub> :0+28+42 Median: 136 <sup>10</sup><br>Standardized IgG against trimeric SARS-CoV-2 spike protein:<br>0+7+14+28+35+42+56 (peak-14 days after the second inoculation) EU: 997.5<br>Multiplex MSD – Anti spike IgG :<br>0+28+42 (peak-14 days after the second inoculation)-(AU/mL): 33830.8<br>Multiplex MSD – RBD IgG:<br>0+28+42 (peak-14 days after the second inoculation)-(AU/mL): 16825.4 | /  | 81.3%  |
|                              |                 |                      | 2. 56-69 years | live SARS-CoV-2 micro-neutralization assay MNA <sub>80</sub> :0+28+42+56 Median: 144 <sup>11</sup><br>Standardized IgG against trimeric SARS-CoV-2 spike protein:<br>0+14+28+35+42+56 (peak-14 days after the second inoculation) EU: 648<br>Multiplex MSD – Anti spike IgG :<br>0+28+42+56 (peak-14 days after the second inoculation)-(AU/mL): 20617   |  |  |
|                              |                 |                      | 3. ≥70 years   | live SARS-CoV-2 micro-neutralization assay MNA <sub>80</sub> :0+28+42+56 Median: 161 <sup>11</sup><br>Standardized IgG against trimeric SARS-CoV-2 spike protein:<br>0+14+28+35+42+56 (peak-14 days after the second inoculation) EU: 525<br>Multiplex MSD – Anti spike IgG :<br>0+28+42+56 (peak-14 days after the second inoculation)-(AU/mL): 19414   |  |  |
| <b>Johnson &amp; Johnson</b> | 1               | /                    | 1. 18-55 years | PRNT:0+28+56+70 (peak-70 days after inoculation) <sup>12</sup><br>(IC50 Log <sub>10</sub> GMT): 321; Seroconversion rate:100%<br>Anti spike IgG:0+28+56+70 (peak-56 days after inoculation)<br>(EU/ml Log <sub>10</sub> GMC): 660; Seroconversion rate:100%  | PRNT-Log <sub>10</sub> GMT:522<br>Anti spike IgG-Log <sub>10</sub> GMC:899   | 72%-United States, 64%-South Africa, 61%-Latin America |
|                              |                 |                      | 2. >65 years   | PRNT:0+14+28 (peak-28 days after inoculation)<br>(IC50 Log <sub>10</sub> GMT): 277; Seroconversion rate:96%<br>Anti spike IgG:0+14+28 (peak-28 days after inoculation)<br>(EU/ml Log <sub>10</sub> GMC): 312; Seroconversion rate:96%  |  |  |



# How the Subunit vaccine Works

| Vaccine developer | Number of doses | Dosing schedule(day) | Trial group | Monitoring time of neutralising antibody(day)   | Human convalescent serum  | Efficacy |
|-------------------|-----------------|----------------------|-------------|---|---|----------|
| <b>Novavax</b>    | 2               | 0+21                 | 18–59 years | PRNT(MN <sub>IC&gt;99%</sub> ):0+21+35 (peak-14 days after the second inoculation) <sup>13</sup><br>GMT:3906<br>Anti-Spike IgG:0+7+21+28+35 (peak-14 days after the second inoculation)<br>(EU/ml): 63,160  | PRNT: 983; S- IgG: 8344<br>Asymptomatic: P-254; S-1661<br>Outpatient symptomatic:<br>P-837; S-7420<br>Hospitalized:<br>P-7457; S-53,391 | 89.30%   |
| <b>ZFSW</b>       | 3               | 0+30+60              | 18–59 years | PRNT <sub>50</sub> :0+44+74 (peak-14 days after the third inoculation) <sup>14</sup><br>GMT: 103; Seroconversion rate: 96.6%<br>RBD specific IgG:0+14+30+44+60+74 (peak-14 days after the third inoculation)<br>GMT: 1782; Seroconversion rate: 99.3% | PRNT <sub>50</sub> - GMT: 51  | Unknown  |





*Do we have any effective methods to detect neutralizing antibodies?*




# Conventional approach- Virus Neutralization Test (VNT)

- “Standard Method”
- High Biosafety requirements
- Time consuming
- Very limited lab resources and high cost



## Determining the PRNT<sub>50</sub>/PRNT<sub>90</sub>

PRNT<sub>50</sub> is the standard endpoint used for PRNT




PRNT<sub>50</sub> is the maximum serum dilution that results in  $\geq 50\%$  reduction in plaque count

Because of these drawbacks of VNT, immunoassays are necessary to vaccine evaluation in both clinical trails and national immunization.

# Immunoassay for detection of neutralizing antibodies

*Neutralizing antibody tests are designed to reflect the protection mechanism in vivo*



Neutralizing antibody is targeting **special sites of S RBD**



S-RBD IgG is targeting **the whole sequence of RBD, BUT only detect IgG type**

Total antibody is target **the whole sequence of RBD, BUT detect IgM, IgG and IgA together.**

*What kinds of serology tests can Mindray provide in the post pandemic era?*




# SARS-COV-2 Antibody Types, Functions and Possible Clinical Scenarios




- The S protein binds to ACE2 through the RBD region of the S1 subunit, mediating **viral attachment** to host cells.
- The FP of SARS-CoV-2 and the two HR domains on S2 are essential for **viral fusion** to host cells.

| Antibody                                      | Main Function   | Possible Clinical Scenario   | Competitors with same design | Progress          |
|---|---|--|------------------------------|-------------------|
| IgG and IgM                                   | IgG and IgM to RBD on S-protein + N protein                     | Aid to diagnosis, indicate disease course  |                              | Already launched  |
| Neutralizing Antibody <i>(Core)</i>           | special sites of S RBD  | May indicate natural infection and vaccine induced immunity  | GenScript                    | Already launched  |
| Neutralizing Antibody <i>(S-RBD IgG)</i>      | the whole sequence of RBD, BUT only detect IgG type             | Indicate past infection<br>May indicate herd immunity<br>May indicate natural infection and vaccine induced immunity | Siemens, Abbott              | Already launched  |
| Neutralizing Antibody <i>(Total antibody)</i> | the whole sequence of RBD, BUT detect IgM, IgG and IgA together | Aid to diagnosis<br>Prevalence screening<br>May indicate natural infection and vaccine induced immunity              | Roche                        | under development |

# IgM & IgG



*Help doctors evaluate the dynamics of the immune response of patients depending on the duration of the disease!*




**Antibody overview and timeline in Sars-CoV-2 infection**

|                    | - Immune response                        |   |  |   |  |  | +                             |  |
|--------------------|--|---|--|---|--|--|-------------------------------|--|
|                    |  |   |  |   |  |  |                               |  |
| <b>Form</b>        | Pentameric                               | Monomeric or Dimeric  |  | Monomeric   |  |  | Monomeric                     |  |
| <b>Heavy chain</b> | ( $\mu$ ) Mu                             | (a) Alpha   |  | (y) Gamma   |  |  | (e) Epsilon                   |  |
| <b>Function</b>    | Humoral immunity<br>Activates complement | Mucosal immunity, Neutralization<br>Fc-mediated effector function |  | Humoral immunity, Neutralization<br>Activates complement, Binds C1q |  |  | Allergy & Parasite<br>defense |  |

*Impact in SARS-CoV-2 unknown*


# Neutralizing antibody tests

The design of NTA assay




competitive method

The design of RBD IgG assay



Indirect method

The design of Total antibody assay




Sandwich method

| Antibody  | Main Function   | Possible Clinical Scenario   | Competitors with same design | Progress          |
|---|---|--|------------------------------|-------------------|
| IgG and IgM                                     | IgG and IgM to RBD on S-protein + N protein                     | Aid to diagnosis, indicate disease course  | Snibe, YHLO                  | Already launched  |
| Neutralizing Antibody ( <i>Core</i> )           | special sites of S RBD  | May indicate natural infection and vaccine induced immunity  | GenScript, Snibe             | launched          |
| Neutralizing Antibody ( <i>S-RBD IgG</i> )      | the whole sequence of RBD, BUT only detect IgG type             | Indicate past infection<br>May indicate herd immunity<br>May indicate natural infection and vaccine induced immunity | Snibe, Siemens, Abbott       | launched          |
| Neutralizing Antibody ( <i>Total antibody</i> ) | the whole sequence of RBD, BUT detect IgM, IgG and IgA together | Aid to diagnosis<br>Prevalence screening<br>May indicate natural infection and vaccine induced immunity              | Roche                        | under development |




# Test Principle-Neutralizing Antibody




competitive method

# Test Principle-S-RBD IgG Assay



Indirect method

# Test Principle-Total Antibody Assay



**Sandwich method**

# Interpretation of Results

1. Non-reactive: S-RBD IgG<10, Total antibody<10, Neutralizing antibody<10
2. Reactive: S-RBD IgG≥10, Total antibody≥10, Neutralizing antibody≥10

As with all analyte determinations, the result should be used in conjunction with information available from clinical evaluation and other diagnostic procedures.

The individual immune response following SARS-CoV-2 infection varies considerably and might give different results with assays from different manufacturers. Results of assays from different manufacturers should not be used interchangeably.




Medicines & Healthcare products  
Regulatory Agency

**WHO International Standard**  
**First WHO International Standard for anti-SARS-CoV-2**  
**immunoglobulin (human)**  
**NIBSC code: 20/136**  
**Instructions for use**  
**(Version 2.0, Dated 17/12/2020)**

## 1. INTENDED USE

The First WHO International Standard for anti-SARS-CoV-2 immunoglobulin is the freeze-dried equivalent of 0.25 mL of pooled plasma obtained from eleven individuals recovered from SARS-CoV-2 infection. The preparation has been evaluated in a WHO International Collaborative study (1). The intended use of the International Standard is for the calibration and harmonisation of serological assays detecting anti-SARS-CoV-2 neutralising antibodies. The preparation can also be used as an internal reference reagent for the harmonisation of binding antibody assays. The preparation has been solvent-detergent treated to minimise the risk of the presence of enveloped viruses (2).




Mindray SARS-CoV-2 S-RBD IgG assay performed good correlation with WHO standard materials.


Pearson correlation coefficient :  $r=0.9978$

Conversion : Mindray AU= $0.8229 \times$ BAU

# WHO Standardization Results of neutralizing antibody and total antibody



NTAb: Pearson correlation coefficient  $r=0.9975$   
 Conversion : Mindray AU= $0.2734 \times$  IU



TAb: Pearson correlation coefficient  $r=0.9995$   
 Conversion : Mindray AU= $1.3062 \times$  BAU

# What can Mindray provide?



- Mindray SARS-CoV-2 **IgM/IgG** assays and
- Mindray SARS-CoV-2 **S RBD IgG** assay
- Intravenous serum or plasma (heparin and citrate), easy to operate;
- Time to first result: **30 minutes**;
- Up to **480 tests/hour** depending on different analyzer models used;
- Fully automatic testing, minimizing infection risk;
- Package: **2 x 50T and 2 x 100T**.

**IgM/IgG** can **support diagnosis** of COVID-19.

**S RBD IgG** can **indicate the immune response induced by natural infection and vaccines**



CL-900i (180T/h)



CL-1000i / CL-1200i(180T/h)



CL-2000i (240T/h)



CL-6000i (480T/h)



# Reference

1. Xia, S. et al. Safety and immunogenicity of an inactivated SARS-CoV-2 vaccine, BBIBP-CorV: a randomised, double-blind, placebo-controlled, phase 1/2 trial. *Lancet Infect Dis* 21, 39-51, doi:10.1016/S1473-3099(20)30831-8 (2021).
2. Xia, S. et al. Effect of an Inactivated Vaccine Against SARS-CoV-2 on Safety and Immunogenicity Outcomes: Interim Analysis of 2 Randomized Clinical Trials. *JAMA* 324, 951-960, doi:10.1001/jama.2020.15543 (2020).
3. Zhang, Y. et al. Safety, tolerability, and immunogenicity of an inactivated SARS-CoV-2 vaccine in healthy adults aged 18-59 years: a randomised, double-blind, placebo-controlled, phase 1/2 clinical trial. *Lancet Infect Dis* 21, 181-192, doi:10.1016/S1473-3099(20)30843-4 (2021).
4. Wu, Z. et al. Safety, tolerability, and immunogenicity of an inactivated SARS-CoV-2 vaccine (CoronaVac) in healthy adults aged 60 years and older: a randomised, double-blind, placebo-controlled, phase 1/2 clinical trial. *Lancet Infect Dis*, doi:10.1016/S1473-3099(20)30987-7 (2021).
5. Jackson, L. A. et al. An mRNA Vaccine against SARS-CoV-2 - Preliminary Report. *N Engl J Med* 383, 1920-1931, doi:10.1056/NEJMoa2022483 (2020).
6. Anderson, E. J. et al. Safety and Immunogenicity of SARS-CoV-2 mRNA-1273 Vaccine in Older Adults. *N Engl J Med* 383, 2427-2438, doi:10.1056/NEJMoa2028436 (2020).
7. Sahin, U. et al. BNT162b2 induces SARS-CoV-2-neutralising antibodies and T cells in humans. medRxiv, 2020.2012.2009.20245175, doi:10.1101/2020.12.09.20245175 (2020).
8. Walsh, E. E. et al. Safety and Immunogenicity of Two RNA-Based Covid-19 Vaccine Candidates. *N Engl J Med* 383, 2439-2450, doi:10.1056/NEJMoa2027906 (2020).
9. Logunov, D. Y. et al. Safety and immunogenicity of an rAd26 and rAd5 vector-based heterologous prime-boost COVID-19 vaccine in two formulations: two open, non-randomised phase 1/2 studies from Russia. *Lancet* 396, 887-897, doi:10.1016/S0140-6736(20)31866-3 (2020).
10. Folegatti, PM. et al. Safety and immunogenicity of the ChAdOx1 nCoV-19 vaccine against SARS-CoV-2: a preliminary report of a phase 1/2, single-blind, randomised controlled trial. *Lancet* 2020; 396: 467–78, doi:10.1016/S0140-6736(20)31604-4 (2020).
11. Ramasamy, M. N. et al. Safety and immunogenicity of ChAdOx1 nCoV-19 vaccine administered in a prime-boost regimen in young and old adults (COV002): a single-blind, randomised, controlled, phase 2/3 trial. *Lancet* 396, 1979-1993, doi:10.1016/S0140-6736(20)32466-1 (2021).
12. Sadoff, J. et al. Interim Results of a Phase 1-2a Trial of Ad26.COV2.S Covid-19 Vaccine. *N Engl J Med*, doi:10.1056/NEJMoa2034201 (2021).
13. Keech, C. et al. Phase 1-2 Trial of a SARS-CoV-2 Recombinant Spike Protein Nanoparticle Vaccine. *N Engl J Med* 383, 2320-2332, doi:10.1056/NEJMoa2026920 (2020).
14. Yang, S. et al. Safety and immunogenicity of a recombinant tandem-repeat dimeric RBD-based protein subunit vaccine (ZF2001) against COVID-19 in adults: two randomised, double-blind, placebo-controlled, phase 1 and 2 trials. *Lancet Infect Dis*, doi:10.1016/S1473-3099(21)00127-4 (2021).